



# InterActive Drill Sequence

## 3.2mmD Straight / Tapered Implant 3.0/3.2



Soft Bone Final



SD2.3  
2.3/2.0mmD

Dense Bone Final



SD2.8  
2.8/2.3mmD

Optional  
Crestal Bone Drills\*



CBD3.2

## 3.7mmD Straight / Tapered Implant 3.0/3.7



SD2.3  
2.3/2.0mmD

Soft Bone Final



SD2.8  
2.8/2.3mmD

Dense Bone Final



SD3.4  
3.4/2.8mmD

Optional  
Crestal Bone Drills\*



CBD3.7

## 4.3mmD Straight / Tapered Implant 3.4/4.3



SD2.3  
2.3/2.0mmD

Optional



SD2.8  
2.8/2.3mmD

Soft Bone Final



SD3.4  
3.4/2.8mmD  
or  
length specific  
NobelReplace®  
RP drill

Dense Bone Final



SD3.8  
3.8/3.4mmD  
or  
length specific  
NobelReplace®  
RP drill

Optional  
Crestal Bone Drills\*



CBD4.2

## 5.0mmD Straight / Tapered Implant 3.4/5.0



SD2.3  
2.3/2.0mmD

Optional



SD2.8  
2.8/2.3mmD



SD3.4  
3.4/2.8mmD

Soft Bone Final



SD3.8  
3.8/3.4mmD  
or  
length specific  
NobelReplace®  
WP drill

Dense Bone Final



SD4.4  
4.4/3.8mmD  
or  
length specific  
NobelReplace®  
WP drill

Optional  
Crestal Bone Drills\*



CBD4.7